



## Pathology Specimen Radiography System

### PathVision from Faxitron®

With the largest field of view and highest resolution detector in the pathology market, PathVision helps you deliver greater efficiency and better patient care.



### PathVision FEATURES

#### High-Resolution Imaging

PathVision has a magnification shelf with field of view guides and automatic position detection that provides up to 10X geometric magnification, which can locate the smallest calcifications in tissue samples to ensure accurate diagnoses.

#### Wide Energy Range

The system can image both thin slices and large specimens, and proprietary algorithms enable optimal image exposure and clarity with Automatic Exposure Control (AEC).

#### Convenient Interface

The continuous surface keyboard allows for ease of use and the system has the ability to save images in a variety of file formats.

#### Easy to Use

No additional training or specialized x-ray requirements are needed to operate PathVision, and the system plugs into any standard A/C outlet.

#### Ergonomic Design

The minimal footprint and elegant form of the PathVision system allows for optimal placement at the point-of-care. Its stainless steel chamber was designed with no right angles to allow for easy cleaning.


#### Improved Efficiency

Routine access to high-resolution imaging with PathVision results in faster final report generation and fewer surgical re-cuts due to a more efficient workflow. The software also has the ability to send multiple annotated images to PACS at the touch of a button.

#### Optimal Imaging Quality

PathVision allows the user to optimally image anything from breast tissue slices and intact mastectomies to bone and fetal remains. The system's large detector has the highest resolution for its size in the industry and accommodates a broad range of samples without sacrificing image quality.

## All Top 10 US Cancer Centers Have Multiple Faxitron Systems



Source Data from US News and World Report - Top Cancer Centers 2017-18. <http://health.usnews.com/best-hospitals/rankings/cancer>

### FAXITRON BIOPTICS, LLC

3440 E Britannia Drive, Suite 150,  
Tucson, Arizona USA 85706

General inquiries: +1.520.399.8180

Toll Free: +1.877.910.0030

Fax: +1.520.399.8182

[sales@faxitron.com](mailto:sales@faxitron.com)

[www.faxitron.com](http://www.faxitron.com)