



Single-Use Wireless Gamma Probe

TruNode® from Faxitron®

TruNode is the first high performance, sterile, single-use, wireless gamma probe that has been optimized for the most common radio-guided surgical procedures. TruNode is cost-effective and delivers superior performance for sentinel lymph node mapping and radioactive seed localization.



TruNode FEATURES

Industry-Leading Resolution

TruNode features a best-in-class 10% energy resolution, enabling the user to detect faint nodes closer to the injection site and facilitating precise hotspot localization.

Fingertip Control

Toggle between three modes of detection (high-sensitivity scan, high-resolution point and I-125 seed) with a simple to use, fingertip control from the sterile field.

Three Detection Modes

With its three detection modes, TruNode optimizes sentinel node, parathyroid and tumor localization procedures.

Minimal Incision Size

With a 10mm gamma probe tip diameter, TruNode is less invasive and easy to maneuver.

Sterile

TruNode is fully sterile, eliminating drapes and reprocessing in the OR.

Voice Reporting

TruNode's compact User Feedback Unit provides both audible and visual user feedback, eliminating the need to monitor a display.

Simple Placement


TruNode has a minimal footprint as the User Feedback Unit clips to any IV pole in the OR.

Single-Patient Use

TruNode was designed to reduce risk of cross-contamination and surgical site infections.

Product Code	Description	UM
120-807605	TruNode Gamma Probe hand piece, single use, sterile	EA
120-807606	TruNode User Feedback Unit console (UFU), reusable Includes: Console, IV Pole Mounting Bracket, Power Charger	EA

All Top 10 US Cancer Centers Have Multiple Faxitron Systems



Source Data from US News and World Report - Top Cancer Centers 2017-18. <http://health.usnews.com/best-hospitals/rankings/cancer>

FAXITRON BIOPTICS, LLC

3440 E Britannia Drive, Suite 150,
Tucson, Arizona USA 85706

General inquiries: +1.520.399.8180

Toll Free: +1.877.910.0030

Fax: +1.520.399.8182

sales@faxitron.com

www.faxitron.com